

Tri-wheel Walkers


Days Aluminium Lightweight Tri-wheel Walker

Light and easy to manoeuvre

Days Aluminium Lightweight Tri-wheel Walkers feature a lightweight aluminium frame that is easy to manoeuvre, making them ideal for use indoors or outdoors.

Key features

- Ergonomic handgrips for optimum comfort
- Adjustable height handles with 5 settings at 25mm intervals
- Lockable loop brakes for security and convenience
- Available with or without a basket
- Includes durable bag to transport and conceal your belongings
- Folds easily for storage and transport using the simple easy to grab pull handle
- ISO 11199-2:2005



A	Handle height 83 - 93cm (32¾ - 36½")
B	Overall depth 61cm (24")
C	Overall width 68.5cm (27")
D	Wheel diameter 19cm (7½")
	Product weight 7.7kg (17 lbs)



Optional basket and tray



*Now ordered in Box Qty

days

With bag only	Description	Qty/ Box
DAYS-TRI-RED-2	Ruby Red	2
DAYS-TRI-BLUE-2	Blue	2
Accessories		
DAYS-091555242	Basket and Tray	

Available colours

