

# IC20 HOMECARE MATTRESS

Middle of the range, the IC20 proves to be popular for most. The medium feel provides excellent support and great bed mobility.

Medical cover's available if required.



## SPECIFICATIONS

<b>Feel Rating</b>	Medium
<b>Mattress Height</b>	200mm
<b>Risk Category</b>	High Risk
<b>User Weight</b>	30kg - 200kg
<b>Warranty</b>	10 Year Warranty



**PRESSURE REDISTRIBUTION**  
The IC20 mattress features the ActiveX™ 84 material is the main comfort layer in this mattress. ActiveX™ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX™ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.



**SOLID SUPPORT FOAM BASE**  
XP Support foam is a high-density foam that is used in the bottom layer of the IC15. XP Support foam is an open cell material that breathes to allow air circulation but is high density to prevent bottoming out. This material will easily contour to any profiling bed.



**TEMPERATURE REGULATING COVER**  
The IC20 cover is made with the Bluetec cooling fabric. This is all way stretch fabric to conform with body and mattress contour. It is under woven with a cooling yarn to regulate body temperature.



**VENTILATED BASE FABRIC**  
Vant fabric is a three-dimensional polypropylene material that has high air permeability and is anti-bacterial. We use this on the side wall of the IC20 mattress covers to greatly increase air movement.



**MEDICAL COVERS AVAILABLE**  
Zerotec fabric is a medical barrier fabric that's has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.

## BED CODES & SIZES

<b>SINGLE</b>	IC2S	1900 X 900mm
<b>LONG SINGLE</b>	IC2LS	2030 x 900mm
<b>KING SINGLE</b>	IC2KS	2030 x 1070mm
<b>DOUBLE</b>	IC2D	1900 x 1350mm
<b>LONG DOUBLE</b>	IC2LD	2030 x 1350mm
<b>QUEEN</b>	IC2Q	2030 x 1520mm
<b>SPLIT QUEEN</b>	IC2SQ	2030 x 760mm (per side)

## CARE & CLEANING

- No flip design reducing manual handling
- Anticrease XP base to allow easy contour for any profiling bed
- Standard cover MVTR: NA. Medical Cover MVTR (g/24h·m2): 1635.0 Average
- Medical Cover BS7175: Crib 5
- Waterlow ≤20
- To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.

## PRESSURE CARE

**Risk Category** High

This rating is based on our average peak pressure of <19mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.