

# Karma Flexx



Only in 16"  
and 18"  
widths

A better  
wheelchair.  
A better fit.



CODE KM-8022



<b>Frame</b>	Light weight aluminium Folding
<b>Footrest</b>	Swing away removable Height adjustable
<b>Armrest</b>	Swing back Height adjustable Ergonomic pads
<b>Wheels</b>	Quick release   24"x 1 3/8" spoke   Solid tyres
<b>Castors</b>	Solid 7"
<b>Brakes</b>	Push to lock
<b>Seat</b>	Height adjustable   Seat extension   Padded
<b>Back</b>	Angle adjustable Padded
<b>Warranty</b>	5 Years - Frame



Robust, stylish and adjustable.

The Flexx is designed with key adjustment to create a suitable sitting environment and to fit individual needs. With ellipse shaped frame design the Flexx is not only stylish but also durable.

The configurable options of the Flexx include adjustments in seat depth, seat height, backrest height, backrest angles, castor, fork and centre of gravity. It can be set up to provide optimum comfort and performance for the user. Heavy duty frame available for 20", 22" and 24".

## OPTIONS

See pages 10-11

## SPECIFICATIONS

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
Flexx 16"x16"	410mm	410mm	500mm	430mm	130kg	14.4kg
Flexx 18"x16"	460mm	410mm	500mm	430mm	130kg	14.6kg
Flexx 18"x18"	460mm	460mm	500mm	430mm	130kg	14.9kg
Flexx 18"x20"	460mm	500mm	520mm	430mm	130kg	15.5kg
Flexx HD 20"x18"	510mm	460mm	520mm	430mm	170kg	17kg
Flexx HD 22"x18"	550mm	460mm	520mm	430mm	170kg	18kg
Flexx HD 24"x18"	605mm	460mm	520mm	430mm	170kg	19kg

